



DIGITAL TOOLS FOR EXTRUSION PROCESS DESIGN

Dr. Josip Matić

Scientific Area Leader, Modeling and Prediction

Pharmaceutical Hot Melt Extrusion



- Hot Melt Extrusion is a one-step, continuous manufacturing step used to produce a variety of pharmaceutical dosage forms, from oral to implantable dosage forms, and filaments for Filament-Based 3D Printing




Formulation perspective

- How to formulate my ASD, extended-release implant, co-extrudate, filament...
- Can my API survive the process?
- Which excipients could be suitable for the extruder and my API?
- How stable will my ASD be?

Process development perspective

- How do I setup the process efficiently for my formulation?
- What quantities of formulation do I need at the early stage?
- How much does my process impact the product quality versus my formulation?
- How to choose the proper process window?
- How to perform scale-up?

HME: Pre-clinical process development




J. Matic, A. Paudel, H. Bauer, R. A. L. Garcia, K. Biedrzycka, and J. G. Khinast, "Developing HME-Based Drug Products Using Emerging Science: a Fast-Track Roadmap from Concept to Clinical Batch," AAPS PharmSciTech, vol. 21, no. 5, p. 176, Jul. 2020.


1D HME




Pharmaceutical Hot Melt Extrusion



Hot Melt Extrusion (HME) is a **continuous, one-step** drug production process.








Pharmaceutical Hot Melt Extrusion





Investigate individual screw elements

Investigate the process and product quality connection

1D HME




18mm → 27mm






Numerous tools for virtual process design



High fidelity simulations

- Simulations of complex fluids, geometries, and free surface flow
- Characterization of the screw elements:
 - What is the conveying capacity and pressure build-up of the screw?
 - How much energy input can be expected?
 - What is the mixing performance of the screw? What kneading or mixing elements should be used?
 - How high is the heat generation?
 - How do elements compare from different vendors and different extruder scales?
- Design and analysis of novel screw elements, tailored to the processing task and formulation







Screw element analysis for HME

SPH screw analysis




**Investigate distributive mixing
with back-pressure**






Types of screw elements


SPH screw analysis




Conveying Elements



Kneading Elements



Mixing Elements






Dimensionless theory

Dimensional Analysis for Extruder Screws (Pawlowski*)

*Pawlowski, J., *Die Ähnlichkeitstheorie in der physikalisch-technischen Forschung*, Berlin/Heidelberg, Springer Publishers, 1971.








$$\frac{\Delta p D}{\eta n L} = A_2 \cdot \left(1 - \frac{1}{A_1} \frac{\dot{V}}{n D^3} \right)$$

For creeping flow ($Re \rightarrow 0$):
linear and
independent of n, D, L, η, ρ
 (geometrical similarity,
 Newtonian, creeping
 flow)

Conveying Elements

Kneading Elements



Mixing Elements





Screw pair performance

Pressure characteristics




How is this useful?

1. Screw configuration assembly and definition of CPPs;
2. Screw transfer during scale-up;




Screw transfer

Screw transfer during scale up





- Define the functional processing section;
- Make sure to transfer those functional processing sections so as to retain their dimensionless length and functionality (conveying capacity, power input, mixing performance...)

18mm  27mm



Screw pair performance


Pressure characteristics




Ongoing Research with Coperion



Behavior of non-Newtonian fluids








Pharmaceutical Hot Melt Extrusion





Investigate individual screw elements

Investigate the process and product quality connection



1D HME



18mm → 27mm




HME processing zones








Numerous tools for virtual process design



1D HME simulations

- Quick and easy-to-use full process simulations
- *In silico* test bench for process development:
 - Test different screw configurations
 - Test different screw speed, throughput, and barrel temperature settings
 - Test different formulations
- Detailed process analysis, e.g. filling degree, melt temperature, residence time distribution...
- *In silico* process setup on different extruders, process transfer, and scale-up

- Input**
- High fidelity simulations;
 - Material data





J. Matić, A. Witschnigg, M. Zagler, S. Eder, and J. Khinast, "A novel in silico scale-up approach for hot melt extrusion processes," Chem. Eng. Sci., vol. 204, pp. 257–269, Aug. 2019.

A. Eitzlmayr et al., "Mechanistic modeling of modular co-rotating twin-screw extruders," Int. J. Pharm., vol. 474, no. 1–2, pp. 157–176, Oct. 2014.

Formulation Properties for Process Setup




Viscosity




Heat Capacity

Guide the Formulation Development to ensure that the Physical Formulation Properties are in line with the Production Technology.

Thermal Conductivity




Spec. Volume






1D HME process design

Coperion ZSK18 process design study



\dot{m} [kg/h]


n [rpm]






1D HME process design

Coperion ZSK18 process design study



\dot{m} [kg/h]

n [rpm]






1D HME process design

Coperion ZSK18 process design study



\dot{m} [kg/h]


n [rpm]






1D HME process design

Coperion ZSK18 process design study



\dot{m} [kg/h]


n [rpm]






1D HME process design

Coperion ZSK18 Process Design Study



\dot{m} [kg/h]



n [rpm]





1D HME process design




Coperion ZSK18 Process Design Study





1D HME process design




Validation







1D HME process design

Validation





Choosing the right process space




18mm

$\dot{m} [kg/h]$ ↓ ↑
 $n [rpm]$ ↓ ↑
 $T_B [^{\circ}C] = const.$
 $\tau [Nm]$

$$SMEC = \frac{2 \cdot \pi \cdot n \cdot \tau}{\dot{m}}$$



Choosing the right process space



18mm

\dot{m} [kg/h]




n [rpm]




T_B [°C] = const.

mRT [s]



Choosing the right process space



18mm

\dot{m} [kg/h]




n [rpm]





T_B [°C] = const.

T_{Die} [°C]

$T_{Barrel} = 150^\circ\text{C}$



Choosing the right process space



18mm

\dot{m} [kg/h]



n [rpm]

T_B [°C] = const.


T_{loc} [°C]

$lmRT$ [s]

$T_{Barrel} = 150^\circ C$



Choosing the right process space




18mm

\dot{m} [kg/h] ↓ ↑

n [rpm] ↑

T_B [°C] = const.

CQA




J. Matić et al., "Towards predicting the product quality in hot-melt extrusion: Small scale extrusion," *Int. J. Pharm. X*, vol. 2, p. 100062, Dec. 2020, doi: 10.1016/j.ijpx.2020.100062.

J. Matić et al., "Towards predicting the product quality in hot-melt extrusion: Pilot plant scale extrusion," *Int. J. Pharm. X*, vol. 3, p. 100084, Dec. 2021, doi: 10.1016/j.ijpx.2021.100084.

Process Scale-up and Transfer


Scale-up Considerations

	10mm	18mm
D [mm]	9.3	18
D_o [mm]	9	17.8
D_i [mm]	5.6	11.8
D_o/D_i [-]	1.6	1.51
Cl [mm]	7.5	15




	27mm
D [mm]	27.2
D_o [mm]	27
D_i [mm]	18.4
D_o/D_i [-]	1.47
Cl [mm]	23



	10mm	18mm	27mm
D_o [mm]	10.5	17.8	27
A_{cr} [mm ²]	65.25	175.32	376.64




Choosing the Right Process Space




$\dot{m} [kg/h]$ ↓ ↑
 $n [rpm]$ ↓ ↑
 $T_B [^{\circ}C] = const.$
 $\tau [Nm]$








Choosing the Right Process Space





$\dot{m} [kg/h]$ ↓ ↑
 $n [rpm]$ ↓ ↑
 $T_B [^{\circ}C] = const.$
 $mRT [s]$



Choosing the Right Process Space



$\dot{m} [kg/h]$ ↓ ↑
 $n [rpm]$ ↓ ↑
 $T_B [^{\circ}C] = const.$
 $T_{loc} [^{\circ}C]$
 $T_{Barrel} = 150^{\circ}C$





HME: *in silico* workflow

Pressure characteristics




Process and material characterization

J. Matic et al., "A novel *in silico* scale-up approach for hot melt extrusion processes," Chem. Eng. Sci., vol. 204, pp. 257–269, Aug. 2019.
J. Matic, A. Paudel, H. Bauer, R. A. L. Garcia, K. Biedrzycka, and J. G. Khinast, "Developing HME-Based Drug Products Using Emerging Science: a Fast-Track Roadmap from Concept to Clinical Batch," AAPS PharmSciTech, vol. 21, no. 5, p. 176, Jul. 2020, doi: 10.1208/s12249-020-01713-0.


Transfer and experimental verification

A


Screw Characterization



Material Properties



B



Virtual process development and scale-up

C




Pharmaceutical HME at RCPE


Holistic product design using ObD





Formulation Design and Product Testing



Process Design & Transfer



In Silico Development Testbed





Pharmaceutical HME at RCPE




Holistic product design using ObD



State-of-the-Art Facilities

Novel Downstream Technologies



Advanced Process Control




Conclusion

J. Matić, A. Paudel, H. Bauer, R. A. L. Garcia, K. Biedrzycka, and J. G. Khinast, "Developing HME-Based Drug Products Using Emerging Science: a Fast-Track Roadmap from Concept to Clinical Batch," AAPS PharmSciTech, vol. 21, no. 5, p. 176, Jul. 2020.

Goal: Ensure high quality medicines developed rapidly, waste free and under the QbD guidelines using *in silico* tools.



Formulation Design




Screw Element Selection

$$\dot{m} [kg/h]$$

$$n [rpm] T_B [^{\circ}C]$$

Processing Design and Transfer





THANK YOU FOR
YOUR ATTENTION!

Dr. Josip Matić

Scientific Area Leader, Modeling and Prediction

